
AVANOS*

**CORGRIP* SR
NG / NI FEEDING TUBE
RETENTION SYSTEM**

Instructions for Use

x 2

x 2

Fig.1

SSL* Clamp Color Reference Guide

8 FR 26-008	EN Pale Blue FR Blue clair DE Hellblau ES Celeste IT Celeste NL Lichtblauw SV Ljusblå	
10 FR 26-010	EN Pink FR Rose DE Rosa ES Rosa IT Rosa NL Roze SV Rosa	
12 FR 26-012	EN White FR Blanc DE Weiß ES Blanco IT Bianco NL Wit SV Vit	
14 FR 26-014	EN Yellow FR Jaune DE Gelb ES Amarillo IT Giallo NL Geel SV Gul	
16 FR 26-016	EN Orange FR Orange DE Orange ES Naranja IT Arancione NL Oranje SV Orange	
18 FR 26-018	EN Tan FR Havane DE Hellbraun ES Tostado IT Marrone chiaro NL Lichtbruin SV Beige	

Fig. 2

 Medical device								Rx Only
		 Choking hazard – Small parts			 Importer		 Unique device identifier	

AVANOS* CORGRIP* SR NG / NI Feeding Tube Retention System

Rx Only: Federal Law (USA) restricts this device to sale by or on the order of a physician.

Description

The AVANOS* CORGRIP* SR NG/NI Feeding Tube Retention System is an accessory to prevent unwanted displacement or removal of NG/NI tubes for adult patients.

⚠ Warning: Not for pediatric use

⚠ Precaution: Use of a securement device may commonly result in minor nosebleed not requiring medical intervention, minor nasal irritation, sinusitis, or patient discomfort. Infrequently occurring factors may also occur as a result of use or misuse of internal securement devices include: abrasion of the nasal mucosa, allergic reaction, contamination, delay in nutrition or medication, nosebleed, nasal perforation, tissue necrosis, aspiration, blockage of airway or nasal passage, or the need for additional medical procedures.

Features

- Catheters: PEBA Material
- Easy to use SSL* (Slot-Slide-Lock) Clamp with French size indication
- Centimeter markings on both catheters for safe insertion
- Biocompatible, long lasting polyurethane monofilament
- Recommended for up to 30 days continuous use

Indications for Use

The AVANOS* CORGRIP* SR NG/NI Feeding Tube Retention System is indicated for use with enteral feeding tubes of 8 FR and greater and NG decompression, suction and drainage tubes up to 18 FR to prevent inadvertent removal or displacement of the tubes for adult patients.

Contraindications for Use

The CORGRIP* SR is contraindicated for patients with facial and/or cranial fractures, nasal airway abnormalities or obstructions. It should also not be used on patients that may pull on their feeding tube or the AVANOS* CORGRIP* SR NG/NI Feeding Tube Retention System to such an extent as to cause serious self-injury.

Should a patient be likely, based on condition or mental status, to continuously pull on the feeding tube or CORGRIP* SR, it is recommended to secure the feeding tube via alternate means to minimize potential injury or tube displacement.

Directions for Insertion

⚠ Do not reuse or reprocess these medical devices. Reuse or reprocessing may 1) adversely affect the known biocompatibility of the device, 2) compromise the structural integrity of the device, 3) lead to the device not performing as intended, or 4) create a risk of contamination and cause the transmission of infectious diseases resulting in a patient injury, illness, or death.

1. It is recommended to perform a nasal exam prior to placement of the CORGRIP* SR to assess for adequate vomer bone support for the CORGRIP* SR Retention System.
2. Position the patient in accordance with facility protocol for feeding tube placement.
3. Remove clear protective sheath from both WHITE and YELLOW* catheters.


Fig. 3

Note: The CORGRIP* SR may be placed before or after the NG/NI tube. For patients with small nares, it may be particularly beneficial to place the CORGRIP* SR prior to tube placement.

4. Thoroughly lubricate the distal tip of both catheters with water soluble lubricant provided in package. Take care to avoid lubricating the magnet tips.
5. Insert the YELLOW* retrieval catheter to the desired depth into desired nostril.

Fig. 4


Use cm markings to guide insertion depth.

Note: Adjust appropriately for smaller patients. Do not exert force when inserting the catheter. If overinserted, damage to sinus cavity may occur.


6. Insert the WHITE catheter into the OPPOSITE nostril to the same distance as the YELLOW* catheter so magnet tips are at the same point behind the vomer bone. Use the printed cm markings on the catheter to match insertion depth of the two catheters.

Fig. 5


7. Some gentle manipulation of the catheters will encourage magnet engagement. When the magnets connect behind the vomer bone, you will hear the click and/or feel the connection. (Reference Progression Summary Figure 2, page 5.)
8. When magnet engagement has occurred, slowly/deliberately withdraw the WHITE catheter until the monofilament is clearly visible and the catheter is completely out of the nostril (approx. 10 cm).

Fig. 6

- Then slowly/deliberately withdraw the YELLOW* catheter until the monofilament is visible and the catheter is completely outside of nostril allowing the advancement of the monofilament behind the vomer bone and out the opposite nostril. (Reference Progression Summary Figures 3–4, page 5.)

Fig. 7

- Now pull the WHITE catheter until the monofilament is completely out of catheter.
- Disconnect the YELLOW* catheter from the monofilament at the magnet. Adjust the monofilament so that an equal length extends from each nostril.
- If the NG/NI tube has not been placed, do so now while taking care not to dislodge the CORGRIP* SR Monofilament.

Securing the CORGRIP* SR to the NG/NI tube with the SSL* (Slot-Slide-Lock) Clamp

Warning: The SSL* Clamp may be a choking hazard when not attached to the feeding tube and/or monofilament.

The SSL* (Slot-Slide-Lock) Clamp allows for manipulation of the clamp well below the patient's mouth. Take care or cover patient's mouth when securing the SSL* Clamp.


Fig. 8

- SLOT the NG/NI tube into the groove of the SSL* Clamp. Then place one or both sides the monofilament onto the flat area of the clamp, in front of the tube.
- SLIDE the clamp up to approximately 1 cm from nostril.
- LOCK the clamp and monofilament into place by pinching the clamp closed.
- Tie a knot with the extra monofilament. See figure 9 step 4 or figure 10 step 2 for details.
Note: Do not tie the knots around the NG/NI tube.
- Trim the excess monofilament about 1 cm below the knot and discard magnet.

Clip both sides of monofilament and tie knot


Fig. 9

OR

Clip one side of monofilament as above and tie loose side below clip.

Fig. 10


Note: French size label on the clamp will be right side up when properly positioned for easy access to open.

Fig. 11

Caution: Take care not to put tension on the monofilament or tube on or close to the nostril as irritation may occur. If irritation occurs remove the CORGRIP SR.

Fig. 12


Note: The SSL* Clamp should not touch the nostril or be secured over the patient's mouth.

Opening the SSL* Clamp

Use the Opening Device provided in the kit if repositioning is desired.

Do not use other objects to open the SSL* Clamp.

Fig. 13

MRI Information

The SSL* Clamp and monofilament are MRI Safe.
The delivery and retrieval catheters are MRI Unsafe.

Caution: Possible interference may occur with pacemaker/cardiac implants due to direct contact with CORGRIP* SR Magnets. Keep the CORGRIP* SR Magnets at least 6 inches/15 cm away from any implanted devices such as pacemakers and defibrillators.

Maintenance

Routinely monitor the position of the SSL* Clamp and tube for indications of migration.

CORGRIP® SR and NG/NI Tube Removal


1. Cut only ONE strand of the monofilament. Gently pull both the CORGRIP® SR and nasal tube out of the nose.
2. If removal of just the CORGRIP® SR is desired, open and remove clamp, cut one end of monofilament and gently remove (tube remains in place).

⚠ **Caution:** Take care to prevent accidental swallowing of the monofilament once the monofilament has been cut.

Reset Monofilament

If connection is not made and procedure needs to be repeated withdraw WHITE catheter. If monofilament is extended beyond the tip, gently pull monofilament from proximal end of WHITE catheter to reset.

Fig. 14


Notification to Users/Patients


Please immediately report any serious incident that may have occurred due to this medical device to the manufacturer or its Authorized Representative and to your National Authority.

Visual of Magnet Engagement Behind Vomer Bone


1. Insert YELLOW* catheter followed by WHITE catheter.


2. Magnets connect.


3. YELLOW* catheter bends to fully engage magnet on tip of each catheter.

4. Monofilament is pulled through WHITE catheter.

5. Gentle removal of catheters exposing monofilament now looped behind vomer bone.

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
Distributed in the USA by Avanos Medical Sales, LLC,
5405 Windward Parkway, Alpharetta, GA 30004 USA.
In USA, 1-844-4AVANOS (1-844-428-2667). www.avanos.com

 Avanos Medical, Inc.,
5405 Windward Parkway,
Alpharetta, GA 30004 USA.

 Avanos Medical Belgium BVBA,
 Leonardo da Vincilaan 1,
1930 Zaventem, Belgium.

Sponsored in Australia by Avanos Medical
Australia Pty Ltd, 475 Victoria Avenue,
Chatswood, NSW 2067, Australia.

COLOR YELLOW* is a Registered Trademark or Trademark of
Avanos Medical, Inc., or its affiliates.

*Registered Trademark or Trademark of Avanos Medical, Inc.,
or its affiliates. © 2018 AVNS. All rights reserved. 2021-09-09
15-M1-1540-02