


Suction Catheters

Atraumatic



Area of application

Suction catheters are used for oro-, nasopharyngeal and tracheobronchial suctioning. Must be attached to a jar/canister via a tube before being connected to vacuum source. The products are used for single use only.

Product range

Description	CH	Connector Color	Tip	Effective Length mm	Packaging	Reference
Suction Catheter, Atraumatic, Straight Tip, with Funnel Connector, 210mm, with 1cm (METRIC) markings						
Atraumatic suction catheter, funnel, metric, straight tip, CH10	10	Black	Straight	210	1/100/600	25210118
Atraumatic suction catheter, funnel, metric, straight tip, CH12	12	White	Straight	210	1/100/600	25210218
Atraumatic suction catheter, funnel, metric, straight tip, CH14	14	Green	Straight	210	1/100/600	25210318
Atraumatic suction catheter, funnel, metric, straight tip, CH16	16	Orange	Straight	210	1/100/600	25210418
Suction Catheter, Atraumatic, Straight Tip, with Funnel Connector, 570mm, with 1cm (METRIC) markings						
Atraumatic suction catheter, funnel, metric, straight tip, CH10	10	Black	Straight	570	1/100/600	25310118
Atraumatic suction catheter, funnel, metric, straight tip, CH12	12	White	Straight	570	1/100/600	25310218
Atraumatic suction catheter, funnel, metric, straight tip, CH14	14	Green	Straight	570	1/100/600	25310318
Atraumatic suction catheter, funnel, metric, straight tip, CH16	16	Orange	Straight	570	1/100/600	25310418
Suction Catheter, Atraumatic, Curved Tip, with Funnel Connector, 570mm, with 1cm (METRIC) markings						
Atraumatic suction catheter, funnel, metric, curved tip, CH12	12	White	Curved	570	1/100/600	25510218
Atraumatic suction catheter, funnel, metric, curved tip, CH14	14	Green	Curved	570	1/100/600	25510318
Suction Catheter, Atraumatic, Straight Tip, with Vacuum Control, 570mm, with 1cm (METRIC) markings						
Atraumatic suction catheter, vacuum control, metric, straight tip, CH10	10	Black	Straight	570	1/100/600	25410118
Atraumatic suction catheter, vacuum control, metric, straight tip, CH12	12	White	Straight	570	1/100/600	25410218
Atraumatic suction catheter, vacuum control, metric, straight tip, CH14	14	Green	Straight	570	1/100/600	25410318
Atraumatic suction catheter, vacuum control, metric, straight tip, CH16	16	Orange	Straight	570	1/100/600	25410418

Sterilization

Sterile, Sterilized by ethylene oxide, Shelf life: 5 years



Suction Catheters

Atraumatic

Certifications

CE-marking	CE 0123
Product classification	Class IIA (EU)
Quality system certified to	ISO 13485, Medical Device Single Audit Program (MDSAP)
European Medical Device Nomenclature (EMDN)	A06010103 CANNULAS AND TIPS FOR APIRATION AND IRRIGATION, SINGLE-USE
Global Medical Device Nomenclature (GMDN)	10749

Instructions for use

Included?	Yes
-----------	-----

Material specifications

Contains latex?	No
Contains phthalates?	No
Contains bisphenols?	No
Contains cotton?	No
Contains paper?	Packaging only
Tube and Connector	Polyvinyl chloride

Packaging information

REF	Primary pack./Single peelpack		Secondary pack./Shelf box		Tertiary pack./Shipper				Pallet Qty. pcs
	Qty. pcs	Dimensions, mm (LxW)	Qty. pcs	Dimensions, mm, (LxWxD)	Qty. pcs	Dimensions mm (LxWxD)	Volume m³	Gross weight kg	
Suction Catheter, Atraumatic, Straight Tip, with Funnel Connector, 210mm, with 1cm (METRIC) markings									
25210118	1	322x33.5	100	325x175x120	600	380x345x370	0.049	7.50	18'000
25210218	1	322x33.5	100	325x175x120	600	380x345x370	0.049	7.50	18'000
25210318	1	322x33.5	100	325x175x120	600	380x345x370	0.049	7.50	18'000
25210418	1	322x33.5	100	325x175x120	600	380x345x370	0.049	7.50	18'000
Suction Catheter, Atraumatic, Straight Tip, with Funnel Connector, 570mm, with 1cm (METRIC) markings									
25310118	1	713x40	100	720x175x120	600	740x375x370	0.103	9.52	9,000
25310218	1	713x40	100	720x175x120	600	740x375x370	0.103	10.54	9,000
25310318	1	713x40	100	720x175x120	600	740x375x370	0.103	11.50	9,000
25310418	1	713x40	100	720x175x120	600	740x375x370	0.103	11.80	9,000
Suction Catheter, Atraumatic, Curved Tip, with Funnel Connector, 570mm, with 1cm (METRIC) markings									
25510218	1	713x40	100	720x175x120	600	740x375x370	0.103	10.54	9,000
25510318		713x40	100	720x175x120	600	740x375x370	0.103	11.50	9,000
Suction Catheter, Atraumatic, Straight Tip, with Vacuum Control, 570mm, with 1cm (METRIC) markings									
25410118	1	713x40	100	720x175x120	600	740x375x370	0.103	10.54	9,000
25410218	1	713x40	100	720x175x120	600	740x375x370	0.103	11.50	9,000
25410318	1	713x40	100	720x175x120	600	740x375x370	0.103	11.80	9,000
25410418	1	713x40	100	720x175x120	600	740x375x370	0.103	10.54	9,000



Suction Catheters

Atraumatic

GTIN Codes			
REF	Primary package	Secondary package	Tertiary package
Suction Catheter, Atraumatic, Straight Tip, with Funnel Connector, 210mm, with 1cm (METRIC) markings			
25210118	6925816600463	16925816600460	26925816600467
25210218	6925816600494	16925816600491	26925816600498
25210318	6925816600524	16925816600521	26925816600528
25210418	6925816600555	16925816600552	26925816600559
Suction Catheter, Atraumatic, Straight Tip, with Funnel Connector, 570mm, with 1cm (METRIC) markings			
25310118	6925816600456	16925816600453	26925816600450
25310218	6925816600487	16925816600484	26925816600481
25310318	6925816600517	16925816600514	26925816600511
25310418	6925816600548	16925816600545	26925816600542
Suction Catheter, Atraumatic, Curved Tip, with Funnel Connector, 570mm, with 1cm (METRIC) markings			
25510218	6925816600579	16925816600576	26925816600573
25510318	6925816600586	16925816600583	26925816600580
Suction Catheter, Atraumatic, Straight Tip, with Vacuum Control, 570mm, with 1cm (METRIC) markings			
25410118	6925816600470	16925816600477	26925816600474
25410218	6925816600500	16925816600507	26925816600504
25410318	6925816600531	16925816600538	26925816600535
25410418	6925816600562	16925816600569	26925816600566

Packaging material	
Primary package	Coated paper, polypropylene/polyethylene film
Secondary package	Single-wall corrugated fiberboard
Tertiary package	Double-wall corrugated fiberboard

Suction Catheters

Atraumatic



Storage and handling	
	THIS WAY UP
	Keep away from rain
	Fragile; handle with care
	Keep away from sunlight
	Do not resterilize
	Do not re-use
	Do not use if package is damaged and consult electronic instructions for use
	Temperature limit (-15° to 49°C)
	Consult instructions for use or use electronic instructions for use

DISTRIBUTED BY

UNOQUIP GMBH
Schützengraben 20
8200 Schaffhausen
Switzerland
info@unoquip.com
www.unoquip.com

SPONSORED IN AUSTRALIA BY

UNOQUIP Pty Ltd
PO Box 1117
Windsor VIC 3181

1300 55 80 79
info@unoquip.au



Well Lead Medical Co.,
Ltd.
C-4 Jinhu Industrial Estate
Hualong, Panyu
Guangzhou 511434
P.R. China



Shanghai Intl. Holding
Corp. GmbH
Eiffestrasse 80
20537 Hamburg
Germany



QS Engineering AG
Rütihof 1
CH 8820 Wädenswil
Switzerland